

Clear Advantage
SPIRIT™
 HYDROCOLLOID ADHESIVE SHEATH

Ordering Information

Size	Order No.	Pack Size
CLEAR ADVANTAGE SPIRIT™ - Style 1		
24mm	5243	30
28mm	5283	30
32mm	5323	30
36mm	5363	30
40mm	5403	30
CLEAR ADVANTAGE SPIRIT™ - Style 2		
24mm	7243	30
28mm	7283	30
32mm	7323	30
36mm	7363	30
40mm	7403	30
CLEAR ADVANTAGE SPIRIT™ - Style 3		
24mm	9243	30
28mm	9283	30
32mm	9323	30
36mm	9363	30
40mm	9403	30

For more information about the CLEAR ADVANTAGE SPIRIT™ Hydrocolloid Adhesive Sheath and other BARD® products, please contact your Bard Representative.



FREEPHONE 0800 0121 699

Fax: 01903 875 085

Email: ukcustomer@crbard.com

Web: www.scripteasy.co.uk



Now part of BARD's Homecare Urology

NURSE ADVICE LINE
FREEPHONE 0800 168 4048

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and directions for use.

Bard Limited, Forest House, Tilgate Forest Business Park, Brighton Road, Crawley, West Sussex, RH11 9BP Tel: +44 (0)1293 527 888 Fax: +44 (0)1293 552 428
 Bard, Clear Advantage, Clear Advantage Spirit, Rochester Medical and Script-easy are trademarks and/or registered trademarks of C. R. Bard, Inc. All other trademarks are the property of their respective owners.
 © 2014 C. R. Bard, Inc. All Rights Reserved. 0514/3579

www.bardmedical.co.uk

BARD | MEDICAL

Clear Advantage
SPIRIT™
 HYDROCOLLOID ADHESIVE SHEATH

Control the level of skin moisture.
 The first and only sheath to use
 hydrocolloid adhesive.



BARD | MEDICAL


Clear Advantage SPIRIT™

HYDROCOLLOID ADHESIVE SHEATH

Dry skin is the key to good adhesion

The hydrocolloid adhesive in CLEAR ADVANTAGE SPIRIT™ silicone allows transpiration across the silicone membrane at a much greater rate - **up to 50% more** than comparable silicone self-adhesive sheaths.

Optimum water vapour transmission




*Data on file at Rochester Medical

Proven to release trapped moisture


Conducts away moisture to **protect skin integrity, minimise skin maceration** and **maintain adhesion**.

Available in a range of styles




Style 1

Standard length sheath with the adhesive positioned centrally. Style 1 sheaths are suitable for men with a good length of penis.



Style 2

Shorter length sheath with the adhesive positioned towards the tip. Style 2 sheaths will fit all sizes of penis and are particularly suitable if the penis is retracted.



Style 3

Standard length sheath with additional adhesive. Style 3 sheaths provide very secure adhesion for men with a good length of penis who are particularly active or have difficulty in keeping other sheaths in place.

