

HydroSil[®]

rose

Female Hydrophilic
Intermittent Catheter

Ordering Information

Size	Order No.	Pack Size
8CH	71408	30
10CH	71410	30
12CH	71412	30
14CH	71414	30
16CH	71416	30
18CH	71418	30

For more information about the HYDROSIL[®] Rose Female Hydrophilic Intermittent Catheter and other BARD[®] products, please contact your Bard Representative.



FREEPHONE 0800 0121 699

Fax: 01903 875 085

Email: ukcustomer@crbard.com

Web: www.scripteasy.co.uk



Now part of BARD[®] Homecare Urology

NURSE ADVICE LINE
FREEPHONE 0800 168 4048

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and directions for use.

Bard Limited, Forest House, Tilgate Forest Business Park, Brighton Road, Crawley, West Sussex, RH11 9BP Tel: +44 (0)1293 527 888 Fax: +44 (0)1293 552 428
Bard, HydroSil, Rochester Medical and Script-easy are trademarks and/or registered trademarks of C. R. Bard, Inc. All other trademarks are the property of their respective owners. © 2014 C. R. Bard, Inc. All Rights Reserved. 0714/3578

www.bardmedical.co.uk

BARD | MEDICAL

HydroSil[®]

rose

Female Hydrophilic
Intermittent Catheter

Soft. Comfortable. Discreet.
The only hydrophilic silicone catheter for ISC.



BARD | MEDICAL



Have you been waiting for a soft and comfortable catheter?

HYDROSIL[®] Rose is made from a **unique hydrophilic silicone** that is soft and flexible yet sufficiently firm and resilient. This makes it the ideal material for a **comfortable 'friction-free'** catheterisation. Easy to fold and so discreet to carry.

Smooth, tapered, seamless tip


Atraumatic insertion

Four comfort sized eyes

Enables optimal drainage

M³ technology

Three layer construction designed to provide optimal handling characteristics while remaining comfortable to use



Pliable innermost layer for easier insertion

Firm middle layer for optimal handling

Ultra-soft outer layer to maximise comfort

Specially designed handle

Provides a firm grip and control over the catheter tip

Ultra-soft hydrophilic silicone

Toxin-free and environmentally friendly. Conforms to the urethral contours and helps towards pain-free catheterisation

Compact size pack

Convenient, discreet and all in one pack

